

The Fume Cupboard Company

Airflow 2 - Educational Ducted Cupboard

The Airflow 2 Educational Ducted Fume Cupboard effectively eliminates all fumes and smoke generated during chemical reactions, providing excellent protection for wet chemistry experiments. It is designed to capture and contain fumes and smoke, and is available in three widths (1000mm, 1200mm, and 1500mm), featuring a chemically resistant Trespa base and a sliding glass sash. This fume cupboard is one of the top performers in the market.

It meets several compliance standards such as Building Bulletin 88 (Rev G9), Energy efficient design and operation of Fume Cupboards Guide 320, and BREAMM Hea 17 and Ene 11. It's worth noting that a suitable extraction system is required for the Fume Cupboard to function correctly. The extract system must meet the requirements of D.W. 154 and Building Bulletin 88.

Furthermore, the Airflow 2 is entirely designed and manufactured in the UK, and it comes in a standard width of 1000mm and a standard height of 1880mm.



If you're interested in Fume Cupboards for Schools, please click the link to find out more.











www.thefumecupboardcompany.co.uk - info@thefumecupboardcompany.co.uk