Airflow 7 Ductless Industrial Fume Hoods.

The Airflow 7, also known as a Ductless Balance Enclosure, is a Class I enclosure that complies with USP 795 requirements for non-sterile environment from hazardous powders and particulates generated on the work surface.

The Purair RX incorporates innovative ductless technology that creates a safe work environment suitable for a wide range of applications, including balance enclosures, compounding, powder mixing, and powder weighing.

Airflow 7 Features & Benefits

- · High efficiency ebm-papst EC blower.
- Energy savings LED lighting.
- Protects the operator from powder and particle hazards.
- Filter clamping eliminates bypass leakage.
- Exhaust canopy allows for thimble ducting to the facility exhaust system.
- Filter blockage alarm.
- Complies with USP 800 USP 795 guidelines.

To learn more about this product, contact us.











